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Model: .....  
 Numer Serii: .....

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Model: .....  
 Serial number: .....

Point of sale stamp

and a legible signature of the Seller: .....

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.  
 03-068 Warszawa, ul. Zawiańska 43  
 Tel. 22 518 36 00, Fax 22 518 36 30  
 www.antar.net  
 e-mail: antar@antar.net

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Date of issue of the manual: April 1, 2022  
 v1-1.04/01/2022

**USER MANUAL WALKER**

**Models: AT02001; AT02002; AT02003**

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**Modelle: AT02001, AT02002, AT02003**

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**HANDLEIDING WALKER**

**Modellen: AT02001; AT02002; AT02003**

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**Modelos: AT02001; AT02002; AT02003**

**MANUALE D'USO CAMMINATORE**

**Modelli: AT02001; AT02002; AT02003**

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**Modeller: AT02001; AT02002; AT02003**

## ENG

Walking sticks are lightweight and easy-to-use products designed to help people with reduced mobility walk. They offer the user support that reduces the risk of falls and allows running with support on both sides. The maximum body weight of the user must not exceed 150 kg.

### MODELS FEATURES:

Depending on their needs and preferences, users can choose from 3 versions of the walkers:

**AT02001** - Version with 4 controllable wheels. The product allows the user to move quickly in all directions and, if necessary, the frame can be locked by turning the brakes. The model is recommended for people with high mobility who do not have balance problems.


**AT02002** – 2-roll version, mutual. The frame can be pushed forward without having to lift it completely.

**AT02003** – Version with 4 adjustable feet, mobile, mutual. This version guarantees maximum stability. For adjustment, the running frame must be lifted completely (in the case of the rigid function) or only on one side (in the mutual function). It is mainly recommended for movements over short distances.

### WARNINGS

- 1) Model AT02001 should not be used by people with balance disorders. This can cause the user to switch with the walking frame.
- 2) Do not lean outside the contour of the walking frame. There is a risk of falls.
- 3) Do not turn with the walker. There is a risk of falls.
- 4) If the 4-wheel frame (AT02001) is not currently needed, your parking brake should always be applied.
- 5) The product must not be used by people weighing more than 150 kg.
- 6) Avoid contact with an open fire. The product is flammable.
- 7) After installation, you need to make sure that the steel bolt on the foot clamp or adjustable rollers is inserted through the holes in both side walls. Otherwise, the product may be unstable and the well may deform. The bolt to change the function from rigid to mutual must also be stuck through the holes in the two side walls.
- 8) The 4-roller model (AT02001) should not be used with a mutual function. It should always be replaced rigidly. To unfold it, hold the hallway frame on the side panels and separate it. After precipitating the frame, the spring-loaded button in the upper central area should fit into the pit and the step frame should be locked when unfolded.

### FOLDING THE WALKING FRAME:

### FOLDING

To fold the balcony, press the latch button and pull the sides of the balcony inwards, and then fold the balcony.

This symbol means: "Maximum user weight"


### APPLICATION

The balcony can be used as a fixed or a crawling option. The rolling option means that the sides of the balcony move in relation to each other allowing the walking frame to be moved forward through lifting only one side. In the fixed option, the walker becomes rigid. It can be moved by lifting the entire walking frame and moving it forward.

### CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

To switch the walking frame from the walking option to the option fixed, remove them from the balcony and install them the corresponding front part. models AT02002 and AT02003 in the constant option cannot be assembled. To him folded, you must switch it to the trailing option.

### WARNING

Before using the product for the first time, adjust its height. To do this, remove the clip, adjust the height - the same for all four legs, and then insert the clip into the appropriate hole. **MAKE SURE THAT THE STEEL PIN PASSES THROUGH THE HOLE ON BOTH SIDES OF THE WALL.** Otherwise, the product may not be stable and the hole may be deformed.

### MOVING WITHIN THE BALCONY

AT02001 - The balcony moves to either side by pushing it. **IF IN THE DATA MOMENT THE USER DOES NOT USE THE PRODUCT, THE BRAKES SHOULD BE TURNED ON.**

AT02002 - The balcony moves forward by lifting its back, pushing it and placing the tail on the ground.

AT02003 - The balcony is moved by lifting it and placing it slightly in front, and in the walking option, by lifting 1 side of the balcony and pushing it forward.

### PRODUCT MAINTENANCE

The product should be periodically cleaned with water and cleaning agents. The product can be treated with disinfectants. When not in use, the product should be stored in a dry room. Avoid prolonged exposure to direct sunlight.

### WARNING:

In case a serious incident that has occurred in relation to the medical device, that directly or indirectly led, might have led or might lead to any of the following:

- a) the death of patient, user or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat;

the above mentioned serious incident should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established;

### WARNING:

In case any pain, allergic reaction or other worrying symptoms, unclear for the user which have occurred in relation to the medical device, the user should consult a healthcare professional;

**NOTE: It is forbidden to use the product contrary to its intended use!**

### SAFE DISPOSAL OF THE MEDICAL DEVICE

When the device has been withdrawn from use it can be treated as municipal waste except for electric medical devices – in such case the electric medical devices shall be disposed in the way described by the law for electrical and electronic waste disposal.

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